



INTERNATIONAL GEMOLOGICAL INSTITUTE

ELECTRONIC COPY

LABORATORY GROWN DIAMOND REPORT

LG508149760

IGI LABORATORY GROWN DIAMOND ID REPORT

01/06/2022

IGI Report Number LG508149760

PEAR BRILLIANT

6.91 X 4.44 X 2.80 MM

Carat Weight	0.51 CARAT
Color Grade	FANCY VIVID BLUE
Clarity Grade	VVS 2
Polish	EXCELLENT
Symmetry	EXCELLENT
Fluorescence	NONE
Inscription(s)	LABGROWN IGI LG508149760

Comments: This Laboratory Grown Diamond was created by Chemical Vapor Deposition (CVD) growth process. Indications of post-growth treatment.

IGI LABORATORY GROWN DIAMOND ID REPORT

01/06/2022

IGI Report Number LG508149760

PEAR BRILLIANT

6.91 X 4.44 X 2.80 MM

Carat Weight	0.51 CARAT
Color Grade	FANCY VIVID BLUE
Clarity Grade	VVS 2
Polish	EXCELLENT
Symmetry	EXCELLENT
Fluorescence	NONE
Inscription(s)	LABGROWN IGI LG508149760

Comments: This Laboratory Grown Diamond was created by Chemical Vapor Deposition (CVD) growth process. Indications of post-growth treatment.

LABORATORY GROWN DIAMOND REPORT

IGI LABORATORY GROWN DIAMOND IDENTIFICATION REPORT

01/06/2022

IGI Report Number LG508149760

Shape and Cutting Style PEAR BRILLIANT

Measurements 6.91 X 4.44 X 2.80 MM

GRADING RESULTS

Carat Weight 0.51 CARAT

Color Grade FANCY VIVID BLUE

Clarity Grade VVS 2

ADDITIONAL GRADING INFORMATION


Polish EXCELLENT

Symmetry EXCELLENT

Fluorescence NONE


Inscription(s) LABGROWN IGI LG508149760

Comments: This Laboratory Grown Diamond was created by Chemical Vapor Deposition (CVD) growth process. Indications of post-growth treatment.



LASERSCRIBESM

Sample Image Used



This Laboratory Grown Diamond (LGD) described in this Report has been analyzed, graded and Laserscribed[®] by International Gemological Institute (IGI). A LGD has essentially the chemical, physical and optical properties as a mined diamond, with the exception of being man-made (a manufactured product). LGD's are typically produced by CVD (chemical vapor deposition) or by HPHT (high pressure high temperature) growth processes and may include post growth modifications to change the color. IGI utilizes the most advanced techniques and equipment currently available including, binocular microscopes, diamond color masters, non-contact optical measuring device, a wide range analytical techniques including FTIR, UV-VIS-NIR, raman spectroscopy, and fluorescence analysis at various excitation wavelengths. This Report includes advanced security features. This Report is neither a guarantee, valuation nor appraisal and by making the report IGI does not agree to purchase or replace the article.

INTERNATIONAL GEMOLOGICAL INSTITUTE, INC.



For Terms & Conditions and to verify this report, please visit www.igi.org