



INTERNATIONAL
GEMOLOGICAL
INSTITUTE

ELECTRONIC COPY

LABORATORY GROWN DIAMOND REPORT

August 4, 2026

IGI Report Number

LG821638345

Description

LABORATORY GROWN DIAMOND

Shape and Cutting Style

PEAR BRILLIANT

Measurements

8.43 X 5.68 X 3.63 MM

GRADING RESULTS

Carat Weight

1.01 CARAT

Color Grade

F

Clarity Grade

VVS 2

ADDITIONAL GRADING INFORMATION

Polish

EXCELLENT

Symmetry

EXCELLENT

Fluorescence


NONE

Inscription(s)

 LG821638345

Comments: As Grown - No indication of post-growth treatment.
This Laboratory Grown Diamond was created by High Pressure High Temperature (HPHT) growth process.
Type II

PROPORTIONS





Sample Image Used

COLOR

D

E

F

G

H

I

J

Faint

Very Light

Light

CLARITY

FL

IF

VVS¹⁻²

VS¹⁻²

SI¹⁻²

I¹⁻³

Flawless

Internally Flawless

Very Very Slightly Included

Very Slightly Included

Slightly Included

Included



© IGI 2020, International Gemological Institute



FD - 10 20

LABORATORY GROWN DIAMOND REPORT



August 4, 2026

IGI Report Number

LG821638345

Description

LABORATORY GROWN DIAMOND

Shape and Cutting Style

PEAR BRILLIANT

Measurements

8.43 X 5.68 X 3.63 MM

GRADING RESULTS

Carat Weight

1.01 CARAT

Color Grade

F

Clarity Grade

VVS 2

ADDITIONAL GRADING INFORMATION

Polish

EXCELLENT

Symmetry

EXCELLENT

Fluorescence


NONE

Inscription(s)

 LG821638345

Comments: As Grown - No indication of post-growth treatment.
This Laboratory Grown Diamond was created by High Pressure High Temperature (HPHT) growth process.
Type II

PROPORTIONS





IGI

August 4, 2026

IGI Report No LG821638345

PEAR BRILLIANT

8.43 X 5.68 X 3.63 MM

Carat Weight

1.01 CARAT

Color Grade

F

Clarity Grade

VVS 2

Depth

63.9%

Table

58%

Graile

Thin To Slightly Thick (Faceted)

Culet

Pointed

Polish

EXCELLENT

Symmetry

EXCELLENT

Fluorescence

NONE

Inscription(s)

 LG821638345

Comments: As Grown - No indication of post-growth treatment.
This Laboratory Grown Diamond was created by High Pressure High Temperature (HPHT) growth process.
Type II

www.igi.org