



ELECTRONIC COPY

LABORATORY GROWN DIAMOND REPORT

August 10, 2026	
IGI Report Number	LG823685339
Description	LABORATORY GROWN DIAMOND
Shape and Cutting Style	ROUND BRILLIANT
Measurements	6.39 - 6.41 X 4.05 MM
GRADING RESULTS	
Carat Weight	1.02 CARAT
Color Grade	F
Clarity Grade	INTERNALLY FLAWLESS
Cut Grade	EXCELLENT


ADDITIONAL GRADING INFORMATION

Polish	EXCELLENT
Symmetry	EXCELLENT
Fluorescence	NONE
Inscription(s)	 LG823685339

Comments: As Grown - No indication of post-growth treatment.
This Laboratory Grown Diamond was created by High Pressure High Temperature (HPHT) growth process.
Type II

LG823685339
Report verification at [igi.org](https://www.igi.org)

PROPORTIONS



Sample Image Used

COLOR

D E F G H I J Faint Very Light Light


CLARITY

FL	IF	VVS ¹⁻²	VS ¹⁻²	SI ¹⁻²	I ¹⁻³
Flawless	Internally Flawless	Very Very Slightly Included	Very Slightly Included	Slightly Included	Included

LABORATORY GROWN DIAMOND REPORT



August 10, 2026	
IGI Report Number	LG823685339
Description	LABORATORY GROWN DIAMOND
Shape and Cutting Style	ROUND BRILLIANT
Measurements	6.39 - 6.41 X 4.05 MM
GRADING RESULTS	
Carat Weight	1.02 CARAT
Color Grade	F
Clarity Grade	INTERNALLY FLAWLESS
Cut Grade	EXCELLENT



ADDITIONAL GRADING INFORMATION

Polish	EXCELLENT
Symmetry	EXCELLENT
Fluorescence	NONE
Inscription(s)	(G) LG823685339
Comments: As Grown - No indication of post-growth treatment.	
This Laboratory Grown Diamond was created by High Pressure High Temperature (HPHT) growth process.	
Type II	



© IGI 2020, International Gemological Institute

FD - 10 20



THIS DOCUMENT WAS PRODUCED WITH THE FOLLOWING SECURITY MEASURES: SPECIAL DOCUMENT PAPER, INK SCREENS, WATERMARK, BACKGROUND DESIGNS, HOLOGRAM AND OTHER SECURITY FEATURES NOT LISTED AND DO EXCEED DOCUMENT SECURITY INDUSTRY GUIDELINES

www.igi.org

August 10, 2026
GI Report No LG823685339
ROUND BRILLIANT

5.39 - 6.41 X 4.05 MM

5.39 - 6.41 X 4.05 MM

Catalytic Weight %	1.02 Catalytic
Color Grade	F

Clarity Grade	I.F.
FL	0.00-0.01
IF	0.01-0.02
VS1	0.02-0.04
VS2	0.04-0.06
VVS1	0.06-0.08
VVS2	0.08-0.10
SI1	0.10-0.12
SI2	0.12-0.14
I1	0.14-0.16
I2	0.16-0.18
I3	0.18-0.20
I4	0.20-0.22
I5	0.22-0.24
I6	0.24-0.26
I7	0.26-0.28
I8	0.28-0.30
I9	0.30-0.32
I10	0.32-0.34
I11	0.34-0.36
I12	0.36-0.38
I13	0.38-0.40
I14	0.40-0.42
I15	0.42-0.44
I16	0.44-0.46
I17	0.46-0.48
I18	0.48-0.50
I19	0.50-0.52
I20	0.52-0.54
I21	0.54-0.56
I22	0.56-0.58
I23	0.58-0.60
I24	0.60-0.62
I25	0.62-0.64
I26	0.64-0.66
I27	0.66-0.68
I28	0.68-0.70
I29	0.70-0.72
I30	0.72-0.74
I31	0.74-0.76
I32	0.76-0.78
I33	0.78-0.80
I34	0.80-0.82
I35	0.82-0.84
I36	0.84-0.86
I37	0.86-0.88
I38	0.88-0.90
I39	0.90-0.92
I40	0.92-0.94
I41	0.94-0.96
I42	0.96-0.98
I43	0.98-1.00
I44	1.00-1.02
I45	1.02-1.04
I46	1.04-1.06
I47	1.06-1.08
I48	1.08-1.10
I49	1.10-1.12
I50	1.12-1.14
I51	1.14-1.16
I52	1.16-1.18
I53	1.18-1.20
I54	1.20-1.22
I55	1.22-1.24
I56	1.24-1.26
I57	1.26-1.28
I58	1.28-1.30
I59	1.30-1.32
I60	1.32-1.34
I61	1.34-1.36
I62	1.36-1.38
I63	1.38-1.40
I64	1.40-1.42
I65	1.42-1.44
I66	1.44-1.46
I67	1.46-1.48
I68	1.48-1.50
I69	1.50-1.52
I70	1.52-1.54
I71	1.54-1.56
I72	1.56-1.58
I73	1.58-1.60
I74	1.60-1.62
I75	1.62-1.64
I76	1.64-1.66
I77	1.66-1.68
I78	1.68-1.70
I79	1.70-1.72
I80	1.72-1.74
I81	1.74-1.76
I82	1.76-1.78
I83	1.78-1.80
I84	1.80-1.82
I85	1.82-1.84
I86	1.84-1.86
I87	1.86-1.88
I88	1.88-1.90
I89	1.90-1.92
I90	1.92-1.94
I91	1.94-1.96
I92	1.96-1.98
I93	1.98-2.00
I94	2.00-2.02
I95	2.02-2.04
I96	2.04-2.06
I97	2.06-2.08
I98	2.08-2.10
I99	2.10-2.12
I100	2.12-2.14
I101	2.14-2.16
I102	2.16-2.18
I103	2.18-2.20
I104	2.20-2.22
I105	2.22-2.24
I106	2.24-2.26
I107	2.26-2.28
I108	2.28-2.30
I109	2.30-2.32
I110	2.32-2.34
I111	2.34-2.36
I112	2.36-2.38
I113	2.38-2.40
I114	2.40-2.42
I115	2.42-2.44
I116	2.44-2.46
I117	2.46-2.48
I118	2.48-2.50
I119	2.50-2.52
I120	2.52-2.54
I121	2.54-2.56
I122	2.56-2.58
I123	2.58-2.60
I124	2.60-2.62
I125	2.62-2.64
I126	2.64-2.66
I127	2.66-2.

Quality Grade	EXCELLENT
Cut Grade	

Depth	63.3%
-------	-------

Table	55%
Girdle	Medium strength

Girdle Medium (Faceted)

Culet	Pointed
-------	---------

	Pointed	EXCELLENT
Culet		
Polish		

Polish	EXCELLENT
Symmetry	EXCELLENT

Fluorescence	NONE
Fluorescence	

Comments: No indication of next growth.

As Grown - No Indication of post-growth treatment.

**This Laboratory Grown Diamond was
Created by High Pressure High**

created by High Pressure High Temperature (HPHT) growth process.